



## Spenco RX® Diabetic Support Foot Beds

Cat. No. 38-141 - Men's

Cat. No. 38-117 - Women's



### Target Consumer:

Consumers with Diabetes who require additional support and cushioning.

### Recommended For:

Part of a total prevention strategy for healthy diabetic feet. Patients should see their doctors and foot care professionals on a regular basis.



### Best For Use In:

Casual and extra depth shoes (or boots).



### Components\*:

EVA, Spenco® Material, Plastazote® Foam,  
Thickness at Heel: 9.58mm (3/8in),  
Thickness at Forefoot: 7.30mm (9/32in)



### Instructions:

Diabetes is a serious illness. The APMA recommends that those who suffer from diabetes should check their feet everyday. Check for areas of irritation (redness), inflammation (swelling), or other changes to their feet. See your healthcare professional at least once a year for a full foot screening and immediately if there are any signs of changes or broken skin or ulcers.

### Features & Benefits:

- 1 **PLASTAZOTE® Foam** provides excellent memory impression and impact absorption to evenly distribute pressure relieving sensitive areas of the foot.
- 2 **Lightweight EVA Base** is body heat activated to mold over time for the precise level of support for the unique shape of your foot. To immediately mold for a full contact custom fit please see your footcare professional. They will use an outside heat source of 230 degrees Fahrenheit.
- 3 **Air Injected Spenco® Material** exclusive cushioning layer provides lasting cushioning, reduced friction and extra comfort.
- 4 **8 Sizes** for correct arch placement for better overall fit, comfort and performance.

**Six Month Unconditional Guarantee**

### Sizing Chart:

MEN'S SIZE CHART				
U.S. mens	9	10	11	12
U.K.	8	9	10	11
European	42	43	44	45

  

WOMEN'S SIZE CHART				
U.S. womens	6	7	8	9
U.K.	4	5	6	7
European	37	38	39	40

**Call toll-free: 1-800-877-3626**  
In Canada: 1-800-387-9538

 [www.spenco.com](http://www.spenco.com)

