




Spenco® Gel Heel Cups

Cat. No. 39-826



Target Consumer:

Consumers seeking contoured cushioning around the heel area.

Consumers suffering from general heel pain, heel spurs or bruises.

Recommended For:

Daily wear – absorbing shock to reduce heel pain.



Best For Use In:

Any dress, casual or athletic shoe.



Components*:

TPR Gel (triple density),

Heel thickness – 9.22mm [3/8"]



Instructions:

Place Spenco® Gel Heel Cup against the inside heel area of shoe. For best fit, unlace shoe and allow maximum room to correctly position foot in the heel cup. Tighten laces after foot is comfortably in place.



Care:

Hand wash with warm water. Air dry.

Features & Benefits:

1

TPR Comfort Gel Thermoplastic Rubber combines superior cushioning and energy return for ultimate comfort.

2

Anatomical Cushioning System
Triple density gels absorb impact and shock where the heel needs it most.

3

Reinforced Heel Cup provides extra cushioning and support for the heel.

One-Year Unconditional Guarantee

Sizing Chart:

Size	Small/Medium	Medium/Large
U.S. Women's	5-8	9-12
U.S. Men's	6-7	8-15
U.K.	4-6	7-14
European	36-40	41-48

Call toll-free: 1-800-877-3626
In Canada: 1-800-387-9538

www.spenco.com



*Measurements are approximate and vary by size.